



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TRANSMED SAHARA III MAXITHERM blood warmer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **SAHARA III** a device for defrosting and heating 6 blood and plasma components packed in plastic bags with a mixing plate,
- **37.0 °C function:**
 - Heats at a constant ambient temperature of 37°C,
 - Heating various blood components,
 - Heated bags with different filling volumes,
- **Characteristics:**
 - Safe way to defrost and heat,
 - There is no risk of contamination with infectious water pathogens often encountered when using a water bath as a defrosting device,
 - Active drying surface of plastic bags provides hygienic environment surrounding blood bags,
 - The temperature of the heated platelet versus the surrounding air is continuously monitored to avoid exceeding the maximum temperature of 37°C of blood components,
 - Standardized defrost and heating method,
 - Slow key response prevents unintentional defrosting of the melt as well as the heating process,
- **MAXITHERM module:**
 - Doubles the load capacity of SAHARA-III to 6 blood bags,
- **Mixing blood components in bags:**
 - Gentle mixing to achieve a homogeneous temperature inside the bag and prevent damage to blood cells,
- **In-built testing module:**
 - Checking instrument functions,
 - Calibration of temperature sensors,
 - No additional measuring devices are required,
- R232 connector,
- English keyboard on the control panel,
- **In the set:**
 - Polish instruction manual (PDF),
 - Power cord,
- Dimensions: 33 x 50 x 33 cm,
- Weight: 14 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

