



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ELEMENT 4 diagnostic-treatment table (New)

- A brand new product,
- Made in Poland,
- ELEMENT 04 stationary rehabilitation tables feature an ABS plastic case, which gives them a modern look and ensures easy maintenance,
- Standard features a three-part recliner with a regulated headrest and adjustment of the top to the Pivot and chair positions, made of plywood,
- Height regulation of the table and the table top to the Pivot and chair positions is performed by an electric actuator,
- The table has a stable base and fully ball-bearing regulation mechanisms,
- There is a free option to order upholstery in the color of your choice, except for "PRESTIGE" colors which are available at extra cost,
- As additional options there are available (at extra cost):
 - Manually regulated armrests,
 - Lowered side shelves,
 - Hanging pad,
 - Mobile wheels system,
 - Foot remote control,
 - Frame control,
- Technical parameters:
 - Length 200 cm,
 - Width 65/70/80 cm,
 - Height: 48 - 100 cm,
 - Headrest angle regulation: -70° to 40°,
 - Backrest angle regulation: 0° to 85°,
 - Tabletop breaking angle regulation: 0° to 30°,
 - Own weight 98 kg,
 - Maximum load 160 kg,
- It has a valid Technical Passport issued,
- 24 months' guarantee for the construction,
- 12 months' guarantee for hydraulic accessories,

