



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS WM-60 endoscopic cart (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- **Real photos of the product,**
- Made in the UK,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1300 VA,
- **OLYMPUS WM-60** was designed to meet the more sophisticated requirements of the endoscopy environment, suitable for transporting endoscopes and medical equipment,
- **Equipped with:**
 - A convenient keyboard holder that can be mounted on either side of the cart,
 - Top dedicated monitor shelf,
 - 3 shelves including 2 regulated in height,
 - Height regulated tripod with 2 hangers,
 - Electrical installation,
 - 4 wheels for necessary mobility including 2 with a brake,
 - Detachable back panel that facilitates cable management,
- **In the set:**
 - Power cords (5 pcs.),
- Dimensions:
 - Overall: 69 x 69 x 127 cm,
 - Monitor shelf: 24 x 24 cm,
 - Upper top: 30 x 57 cm,
 - 3 shelves: 54 x 50 cm,
- Maximum load: 140 kg,
- Weight: 75 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

