



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **OLYMPUS WM-60 endoscopic cart (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- **Real photos of the product,**
- Made in the UK,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1300 VA,
- **OLYMPUS WM-60** was designed to meet the more sophisticated requirements of the endoscopy environment, suitable for transporting endoscopes and medical equipment,
- **Equipped with:**
  - A convenient keyboard holder that can be mounted on either side of the cart,
  - Top dedicated monitor shelf,
  - 3 shelves including 2 regulated in height,
  - Height regulated tripod with 2 hangers,
  - Electrical installation,
  - 4 wheels for necessary mobility including 2 with a brake,
  - Detachable back panel that facilitates cable management,
- **In the set:**
  - Power cords (5 pcs.),
- Dimensions:
  - Overall: 69 x 69 x 127 cm,
  - Monitor shelf: 24 x 24 cm,
  - Upper top: 30 x 57 cm,
  - 3 shelves: 54 x 50 cm,
- Maximum load: 140 kg,
- Weight: 75 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

