



OOSTWOUD diagnostic - treatment table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the Netherlands,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- 2 - segment diagnostic - treatment Oostwoud table,
- **Medically approved upholstery in black after purchase will be replaced with new one to the color of customer's choice (from our color pattern book),**
- **Specification:**
 - Frame made of powder-coated profiles,
 - Electric height regulation by a wired remote control in the range: 54-90 cm,
 - Manual headrest angle regulation in the range: 0 - 55°,
 - Paper pad hanger 55 cm wide,
 - Mobile on 15 cm diameter wheels with a central brake and directional travel option,
- **In the set:**
 - Power cord,
 - A remote control,
- **Dimensions:**
 - Bedspace: 195 x 75 cm,
 - Overall: 195 x 75 x 54/90 cm,
- Weight: 85 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

