



## **MONTANA KETTLER stepper (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany (2017),
- Power supply via 2 AA type batteries,
- **Montana stepper by Kettler** will work well for both amateur workouts and intense exercises of people who are a bit more advanced in terms of training. The device features load adjustment based on 12 levels,
- **Workout computer functions:**
  - Time measurement,
  - Measuring the number of steps covered in a training unit,
  - Measurement of pace,
  - Total altitude measurement,
  - Indoor temperature measurement,
  - Energy consumption/heart rate measurement,
  - Resting heart rate measurement with fitness assessment,
- **English menu,**
- **German keyboard on the control panel,**
- **In the set:**
  - Kettler Cardio Heart Rate Clip,
- Dimensions: 95 x 80 x 155 cm,
- Maximum load: 110 kg,
- Weight: 30 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

