



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EMS MEDI-LINK Model 70 shortwave diathermy (electrosurgical unit) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in England,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 200 VA,
- It has an ultrasonic, interference and laser module,
- Frequency of emitting waves: F= 27.12 MHz,
- Peak power: 0 to 50 W,
- Pulse repetition rate: 5 to 400 Hz,
- Pulse width: 20 to 400 µs,
- Treatment duration: 0 to 30 min,
- Built-in display: 240 x 128 pixels,
- Foiled touchscreen control panel,
- Entirely placed on a mobile table:
 - The table has a pull-out top with dimensions: 23 x 26,
 - A drawer with dimensions: 19 x 24.5 x 13 cm,
 - Mesh basket with dimensions: 54 x 29 x 10 cm,
- **Menu to be chosen in: EN, DE, FR,**
- **English keyboard on the control panel,**
- **In the set:**
 - Head (electrode) with dimensions: diameter 17 cm, length 12 cm, 10W,
 - Laser head - wavelength 850 nm, output power 100 mW,
 - Laser head - wavelength 640 nm, output power 10 mW,
 - LED head - wavelength 950 and 640 nm, output power 10 x 15 mW and 9 x 10 mW,
 - 2 and 4 pole cable (no electrodes included),
 - Protective goggles (2 pcs.)
 - Power cord,
- Dimensions: 67 x 45 x 143 cm,
- Weight: 28 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

