



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Endoscopic cart (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **Equipped with:**
 - Swivel top (360°) for monitor with dimensions: 45 x 40 cm,
 - Pull-out keyboard shelf with dimensions: 35 x 39 cm,
 - 4 shelves with dimensions: 44 x 45 cm,
 - A drawer with dimensions: 35 x 34 cm,
 - Side panel for power cords,
 - Power switch on the side panel,
- Made of powder coated steel,
- Mobile on 4 wheels with a diameter of 12 cm including 2 with a brake,
- **In the set:**
 - Power cords (4 pcs.),
- Overall dimensions: 56 x 61 x 147 cm,
- Weight: 55 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

