



EMED ES 300 diathermy (electrosurgical unit) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 600 VA,
- **EMED ES300** is a universal diathermy operating on the basis of high-frequency current flow. It allows performing a variety of procedures in general surgery, vascular surgery, gynecology, oncology and others.
- **It has the following functions:**
 - Clear scission - 250 W dla 300 Ω,
 - Coagulation scission- 150 W for 300 Ω,
 - Soft coagulation - 150 W for 150 Ω,
 - Forced coagulation - 150 W for 150 Ω,
 - Bipolar coagulation - 99 W for 40 Ω,
- Smooth power regulation,
- Digital displays of set power,
- Linear display - control of neutral electrode application,
- Foiled control panel,
- **In the set:**
 - Instruction manual in: PL, EN (PDF),
 - A foot switch,
 - Brand new cable with BOWA clip for neutral electrodes - **no electrodes included, can be added at extra cost,**
 - BOWA monopolar electrode holder (two buttons, 2.4 mm shaft), **(no electrodes included)**
 - Brand-new BOWA bipolar electrode lead **(can be equipped with bipolar tweezers at additional cost),**
 - Power cord,
- Dimensions: 34 x 37 x 16 cm,
- Weight: 8 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

