



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK RSL 2100 Slit lamp - white (New)

- New,
- Made in Germany,
- Advantages and features of Rodenstock RSL 2100 Slit Lamp:
 - Zeiss type lamp (bottom-light),
 - LED illumination,
 - Stereomicroscope with five magnification ranges,
 - High quality optics,
 - 5-step magnification change,
 - Rich set of filters and accessories,
 - All standard filters: green, cobalt, gray,
 - Complete with chinrest, transformer, accessories,
- The slit lamp is equipped with a complete set of filters, such as:
 - Heat absorbing filter,
 - ND filter,
 - Red free filter ,
 - Cobalt blue filter,
- Additional options (at an extra charge);
 - Yellow filter,
 - Applanation tonometer,
- Technical specifications:
 - Microscope:
 - Type: Galilean,
 - Version with five magnifications,
 - 12.5x magnification eyepiece,
 - Total magnification factor: 6x, 10x, 16x, 25x, 40x,
 - Pupil spacing: 52mm~80mm,
 - Diopter adjustment: -8 D to +8 D,
 - Field of view: 36.2, 22.3, 14, 8.9, 5.7 mm,
 - Slit illumination:
 - Slit width: Infinitely adjustable 1-14 mm,
 - Slot length: Infinitely adjustable 0-14 mm,
 - Slot diameters: Ø14 mm, Ø8 mm, Ø3.5 mm, Ø0.2 mm,
 - Slot angle adjustment: 0-180 degrees,
 - Filters: yellow, red-free, neutral, blue, IR heat-absorbing filter, diffuser,
 - Luminance: ≥150klx,
 - Chin:
 - Vertical movement: 80mm,
 - Dimension and weight:
 - Dimension (WDH): 330 x 390 - 460 x 460 mm,
 - Gross weight: 16 kg,
 - Power requirements:
 - Power input: 100 - 240 V, 50/60 Hz,
 - Power output: 12V - 15V / 4A,
- It has a current Technical Passport issued,
- Warranty 12 months, Valid technical passport issued,

