



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **RODENSTOCK RSL 2600 Slit lamp - white (New)**

- New,
- Made in Germany,
- Advantages and features of Rodenstock RSL 2600 Slit Lamp:
  - Zeiss type lamp (bottom-lighting),
  - LED illumination,
  - Stereomicroscope with five magnification ranges,
  - High quality optics,
  - 5-step magnification change,
  - Rich set of filters and accessories,
  - All standard filters: green, cobalt, gray,
  - Complete with chinrest, transformer, accessories,
- The slit lamp is equipped with a complete set of filters, such as:
  - Heat absorbing filter,
  - ND filter,
  - Red free filter ,
  - Cobalt blue filter,
- Additional options (at an extra charge):
  - Yellow filter,
  - Applanation tonometer,
- Technical specifications:
  - Microscope:
    - Type: Galilean,
    - Version with five magnifications,
    - 12.5x magnification eyepiece,
    - Total magnification factor: 6x, 10x, 16x, 25x, 40x,
    - Pupil spacing: 52mm~80mm,
    - Diopter adjustment: -8 D to +8 D,
    - Field of view: 36.2, 22.3, 14, 8.9, 5.7 mm,
  - Slit illumination:
    - Slit width: Infinitely adjustable 1-14 mm,
    - Slot length: Infinitely adjustable 0-14 mm,
    - Slot diameters: Ø14 mm, Ø8 mm, Ø3.5 mm, Ø0.2 mm,
    - Slot angle adjustment: 0-180 degrees,
    - Filters: yellow, red-free, neutral, blue, IR heat-absorbing filter, diffuser,
    - Luminance: ≥150klx,
  - Chin:
    - Vertical movement: 80mm,
- Dimension and weight:
  - Dimension (WDH): 330 x 390 - 460 x 650 mm,
  - Gross weight: 17 kg,
- Power requirements:
  - Power input: 100 - 240 V, 50/60 Hz,
  - Power output: 12V - 15V / 4A,
- It has a current Technical Passport issued,
- Warranty 12 months,

