



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **RODENSTOCK PHOROVIST 800 Phoropter (New)**

- New,
- Made in Germany,
- Advantages and features of the PHOROVIST 800 Phoropter:
  - Simple, intuitive operation,
  - Very robust design,
  - Fast and precise,
  - Lenses with multilayer coating,
  - Auxiliary lenses,
  - Cross cylinder + rotary prism,
  - Convergence system,
- Specifications:
  - Measurement range:
    - Sphere:
      - -19.00 to +16.75 D (steps in 0.12, 0.25 D increments),
    - Cylinder:
      - 0.00 to -6.00 D (steps in 0.12, 0.25 D increments),
      - 0.00 to -8.00 D (when optional lens is used)
    - PD: 50 to 80 mm,
    - Prism: 0 to 20  $\Delta$  (steps of 1  $\Delta$ ),
    - Cylinder axis 0 to 180  $^{\circ}$ ,
  - REFRACTOR HEAD:
    - Auxiliary lens: Occluder, pinhole, red/green filter, polarizing filter (45 $^{\circ}$ /135 $^{\circ}$ ), fixed cross-cylinder lens ( $\pm 0.50$  D), red/white Maddox rod vertical/horizontal, prism (6, 10  $\Delta$  I, U), +0.12 D lens,
    - Crossover cylinder:  $\pm 0.25$  D,
    - Additional lenses: -0.12 D cylindrical lens, -2.00 D,
    - Retinoscopy: +2.00 D (50 cm),
    - Effective field of view: 18.5 mm,
  - Dimensions and weight:
    - WDH dimensions: 340 x 110 x 310 mm,
    - Weight: Approx. 4.8 kg,
- Has a current Technical Passport issued,
- Warranty 12 months,

