



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK CV 800 P LCD panel (New)

- New,
- Made in Germany,
- Advantages and features of the RODENSTOCK CV 800 P LCD panel:
 - Wi-Fi/Bluetooth connection with Phoromat 2000,
 - Accurate linear polarization,
 - Polarization tests including heterophoria,
 - Red-green optotypes for binocular vision tests,
 - Color vision tests,
 - Modern solution for ophthalmology offices,
 - Convenience of use,
 - Wall mounting,
- Specifications:
 - Display: 24" Full-HD LCD panel,
 - Distance units: feet or meters,
 - Working distance: 2 to 7 m (measuring step 0.1 m),
 - Program: 2 programs with up to 20 optotypes,
 - External input/output port: 1 x USB 2.0, RS-232C,
 - Communication: Bluetooth/WiFi for wireless connection, RS-232 for wired connection,
 - Remote control: included,
 - AVAILABLE OPTOTYPES AND TESTS:
 - Standard: Snellen, Landolt, letters, numbers, characters,
 - Color vision: Color vision tables,
 - Red - Green: Schober, Worth, Bichrom-Balance, Fixation Disparity,
 - Low vision: ETDRS, Crowding Bars, Amsler,
 - Contrast: Contrast Adjustment,
 - Children: HOTV, hands, images, symbols
 - Additional: Patient education,
 - DIMENSIONS AND ELECTRICAL REQUIREMENTS:
 - WDH dimensions: 583 x 55 x 423 mm,
 - Weight: Approx. 5.7 kg,
 - Power supply: AC 100 V to 240 V,
 - Frequency: 50/60 Hz,,
 - Power consumption: Maximum 35 - 45 W,
- Has a current Technical Passport issued,
- Warranty 12 months,

