



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TopaScope Non-contact tonometer with pachymeter (New)

- New,
- Made in Germany,
- Product advantages and features:
 - Automatic 3D tracking and automatic one-touch measurement,
 - IOP measurement using soft air blast without contact,
 - CCT measurement using Scheimpflug technology,
 - Pachymetry for IOP correction,
 - Excellent repeatability of measurements,
 - Internal patient database,
 - Very fast and precise measurement,
 - Automatic positioning on the pupil,
 - Very user-friendly interface on 10.1" touch panel,
 - Built-in thermal printer with auto-cut function,
 - Advanced, state-of-the-art technology,
- Technical data:
 - INTRAOCULAR PRESSURE MEASUREMENT:
 - Measurement range: 1 - 60 mmHg,
 - Measurement range setting: Auto / 30 mmHg / 60 mmHg,
 - Measurement principle: Blowing method,
 - Displayed units: mmHg / hPa,
 - Working distance: 11 mm,
 - PACHYMETRIC MEASUREMENT:
 - Measurement range: 400 - 800 μ m,
 - Measurement principle: Slit image on the central cornea,
 - Light source: blue LED,
 - MAIN UNIT:
 - Measurement modes: Fully automatic / automatic / manual,
 - Alignment: Fully automatic 3D tracking,
 - Chin: automatic,
 - Display: 10.1-inch LCD touch screen,
 - Printer: Thermal linear printer,
 - Range of motion:
 - Front / Back: 40 mm,
 - Left / Right: 90 mm,
 - Top / Bottom: 30 mm,
 - Chin movement range: Up / Down: 70 mm,
 - Interface: USB, RS-232,
- DIMENSIONS AND ELECTRICAL REQUIREMENTS:
 - Dimensions: 282 mm (W) × 500 mm (D) × 500 mm (H),
 - Weight: 17 kg,
 - Power supply: AC100 V to 240 V (1.1 A~1.9 A) 50/60 Hz,
- Has a current Technical Passport issued,
- Warranty 12 months, Valid technical passport issued,

