



RODENSTOCK R30 Nyctometer (New)

- New,
- Made in Germany,
- Product advantages and features:
 - Simple and fast operation,
 - Small and compact housing,
 - Modern and precise technological solutions,
 - Test set for occupational medicine,
 - Comprehensive testing for drivers,
 - Innovative solutions,
 - Visual acuity tests for distance and near vision,
 - Standardized tests for drivers and Ishichara,
 - Pediatric test with single optotype,
 - Matsubara test,
 - Twilight vision tests with and without glare,
 - Portable and easy to use,
 - Easy to operate with a 7-inch tablet,
 - Carrying bag,
- Specifications:
 - TESTS:
 - Test distance: 33 cm, 40 cm, 57 cm, 100 cm, unlimited,
 - Test field size: 5°,
 - Photopic test field brightness: 120 cd/m,
 - Mesopic test field brightness without glare: 0.032 cd/m,
 - Mesopic field brightness test with glare: 0.1 cd/m,
 - PROTECTION:
 - Protection level: EMVG / DIN EN 60601-1,
 - Protection class: MPG 1,
 - DIMENSIONS AND ELECTRICAL REQUIREMENTS:
 - WDH dimensions: approx. 265 x 330 x 460 mm,
 - Weight: Approx. 9.8 kg,
 - Voltage: 230 V,
 - Frequency: 50/60 Hz,
 - Rated input voltage: 230 V, 50/60 Hz,
 - Power consumption: 40 VA,
 - Glare intensity: 0.35 lx,
- Has a current Technical Passport issued,
- Warranty 12 months,

