



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **RA LAMB E66/3 Laboratory cooling plate (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in England,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 450 W,
- LAMB E66 is a simple and efficient laboratory cooling plate, ideal for cooling paraffin blocks,
- Dimensions of the plate: 37 x 29.5 cm,
- Digital display,
- Cooling temperature control range from 0° to -10°,
- Has a pull-out drawer with dimensions: 37 x 30 x 3 cm,
- Includes:
  - Power cord,
- Overall dimensions: 38 x 61 x 25 cm,
- Weight: 20 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,
- Valid technical passport issued,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

