



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**PHILIPS HD15 Ultrasound scanner (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in USA, (2012),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 750 VA,
- PHILIPS HD15 ultrasound scanner features excellent imaging quality and high performance. With advantages useful for most everyday examinations and a high level of diagnostic confidence to meet the toughest challenges, the word "value" takes on a new meaning with this instrument. The HD15 is designed to streamline and accelerate the performance of diagnostically confidence-enhancing ultrasound examinations that are a source of clear data, especially for patients and pathologies that are "technically difficult to image."
- Imaging modes:
  - 2D, 3D,
  - M,
  - Color Doppler,
  - Power Doppler,
  - Pulse Wave (PW),
  - Dual mode,
  - CW,
  - Harmonic imaging,
  - Zoom,
- Key advantages:
  - High-resolution A/D technology, 4x parallel processing and Microfine EX focusing,
  - SONOCT real-time complex imaging,
  - Advanced adaptive XRES image processing,
  - Highly sensitive contrast-enhanced ultrasound,
  - Clearer pure wave crystal technology for imaging a wider range of patients,
  - Tissue-specific imaging presets and iSCAN button optimization in 2D, Doppler and color,
  - Active native data allows manipulation of examination parameters and image settings
  - Tissue-specific imaging presets and iSCAN button optimization in 2D, Doppler and color,
  - Active native data allows manipulation of examination parameters and image settings,
- Equipped with:
  - HD15 Revision 3.0.3 system,
  - Monitor tilting 20" LCD screen,
  - Backlit control panel with built-in alphanumeric keypad and height adjustment,
  - 3 head ports,
  - ECG port,
  - USB ports,
  - Mobile on four casters including two with brake,
- Menu in language: English,
- Keypad on control panel in language: English,
- Includes:
  - Sony UP-D897 Videoprinter,
  - ECG 3-lead cable,
  - Linear head L12-3,
  - Sector head S5-2,
  - Power cord,
- Dimensions: 57 x 110 x 150 cm,
- Weight: 100 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN DUTCH  
Мы говорим по русски  
Мы говорим по украински

