



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**OHMEDA GIRAFFE Neonatal incubator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- OHMEDA GIRAFFE neonatal incubator is an incubator that promotes the growth of fragile newborns and allows for intensive neonatal care - thanks to its high thermal efficiency, innovative design and ease of adjusting temperature, humidity, lighting, noise level or oxygen content, it is possible to stabilize the health of patients,
- Equipped with a medical approved mattress in cream color,
- It has built-in systems:
  - Maintaining constant temperature,
  - Regulation of air humidity,
  - Alarm,
- Equipped with:
  - Digital displays for current and desired skin and ambient temperatures,
  - 5 nursing openings for easy access to the patient,
  - Distilled water container,
  - Pull-out tray,
  - Pull-out drawer on both sides,
  - Grommet seals for wires and drains - 13 pieces,
  - Two sockets for connecting two skin temperature sensors (one sensor included),
  - Foiled control panel,
  - Electronic scale,
  - Power sockets for connecting additional devices - 3 pcs,
- Specifications:
  - Infinitely variable manual tilt of the recliner to Trendelenburg and anti-Trendelenburg positions within 12º,
  - Electric height adjustment of the incubator by foot switch in the range: 85 - 112 cm,
  - Adjustable humidity of the air in the chamber in the range: 30 - 95 %,
  - Air temperature: 20 - 37º C,
  - Skin temperature: 32 - 37 º C,
  - Oxygen: 21 - 65%,
  - Heating time: < 35 min,
  - Fresh air flow rate: 30 L / min,
- Mobile on 4 wheels with Ø 12 cm, all with brakes,
- Menu in language: English,
- Included:
  - User manual in English (PDF),
  - Skin temperature sensor,
  - Air hose,
  - Power cord,
- Dimensions:
  - General: 110 x 65 x 147/173 cm,
  - Chamber: 82 x 54 x 35 cm,
  - Mattress: 48 x 65 cm,
  - Drawer: 52 x 48 x 20 cm,
- Weight: 140 kg,
- It has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN DUTCH  
Мы говорим по русски  
Мы говорим по украински

