



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FOSHAN ZC-S300 Dental chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in PRC,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 800 VA
- Upholstery with medical certificate in blue color - it has signs of use (it will be replaced with a new one, in the color chosen by the customer, from our color chart),
- Specifications:
 - Rigid frame and headrest construction,
 - Ergonomic 3D seamless upholstery,
 - No need to bolt the seat to the floor,
 - Armrest can be mounted on the right or left side of the seat,
 - Infinite electric adjustment by foot switch located in the base of the chair:
 - Seat height in the range: 50 - 85 cm,
 - Backrest angle in the range of 0 - 80°,
 - Manual adjustment of headrest height,
 - Optimized to work with multiple units,
 - Equipped with an arm-mounted lamp with manual adjustment in multiple planes,
- Dimensions: 140 x 82 x 190 cm,
- Weight: 145 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

