



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EVA Radiofrequency laser (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Italy, (2017),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 100 W,
- EVA radiofrequency laser is the first device in the world that uses a unique advanced VDM (Vaginal Dynamic RadioFrequency) technology - it is based on Quadripolar Dynamic Energy, which guarantees high accuracy of treatments, on specific layers of vaginal tissue. Thanks to the system of controlled energy emission, the risk of overheating is reduced, eliminating the risk of burning,
- The RSS (RadioFrequency Safety System) allows the full potential of VDR technology to be unleashed, guaranteeing safety during each session. The system includes: handpiece movement monitoring system, 4-point temperature monitoring system, bio-feedback button, allowing the patient to interrupt the treatment,
- UPR - UltraPulsed Radioporation system for transdermal delivery of active substances (e.g. vitamins, hyaluronic acid, hormones) to vaginal and vulvar tissue,
- Automatically controlled electrodes and monitoring of the handpiece's motion sensors prevent the therapist from staying in one place too long, protecting the tissues from overheating,
- In addition, the temperature is constantly monitored to make sure that the therapeutic level is reached on the one hand, and that the upper limit is not exceeded on the other,
- Advantages:
 - Effective endovaginal and labial treatments,
 - A lunchtime procedure that lasts 20 minutes,
 - Non-invasive walk-in/ walk-out procedure,
 - Continuous monitoring of all key treatment parameters,
 - Medical certification,
- Indications:
 - Lack of or reduced sexual satisfaction,
 - Vaginal and vulvar atrophy (VVA),
 - Symptoms of Menopausal Genitourinary Syndrome (GSM),
 - Vaginal relaxation syndrome,
 - Mild stress urinary incontinence,
 - Willingness to rejuvenate the vagina and labia,
- Specifications:
 - RF power: 55W ± 10%,
 - Power control: digital PWM 5% - 100%,
 - VDR technology: 1 Mhz quad dynamic RADIO4,
 - RSS Technology: RF multiple safety control technology,
 - Programs: Software with pre-saved programs,
 - Controls: Rotary encoder with parameter selection button, touch screen, foot switch, Bio-feedback,
 - Security: Power switch with key, emergency button, Bio-feedback, RSS technology,
 - Connectivity: LAN, Bluetooth, USB,
 - Mobile on cart with 3 wheels, all with brake,
 - Selectable menu in language: EN, IT,
- Included:
 - User manual in Polish,
 - Head with cable + caps (6 pcs - 3 internal and 3 external).
 - Foot switch,
 - Bio-feedback button,
 - Power cord,
- Dimensions: 45 x 54 x 113 cm,
- Weight: 18 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3 months for the international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

