



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**DRAGER ISOLETTE 8000 Neonatal incubator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in USA, (2012),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Isolette 8000 allows continuous monitoring of central and peripheral body temperature of newborns, is an ideal incubator for neonatal intensive care units,
- Equipped with a new medical approved mattress in cream color,
- Advantages of Isolette 8000:
  - Dual air curtains and an advanced servo-regulation system maintain stable conditions inside the incubator during patient access,
  - The thermomonitoring system gives the operator the ability to observe the baby's central and peripheral skin temperature, also in servo control mode,
  - Trends are clearly presented on the LCD display,
  - Foiled control panel,
  - Easy access to the newborn,
  - Convenient electronic height adjustment of the incubator,
  - Manual adjustment with a knob of the chamber angle (Trendelenburg or anti-Trendelenburg position) within the range:  $\pm 12^\circ$ ,
  - Low weight, compact design,
- Equipped with:
  - Tilting front and rear access panel,
  - 6 nursing openings for easy access to the patient,
  - Grommet seals for wires and drains,
  - Socket for connecting two skin temperature sensors,
  - Rails for attaching additional accessories,
  - X-ray cassette drawer under the incubator chamber,
  - 3 storage drawers,
  - Self-adjusting casters, and ergonomic foot brakes allow the incubator to be moved easily, even in limited ward space,
  - RS-232 output,
  - Humidifier with a capacity of 1600 ml,
  - Bottle holder,
- Specifications:
  - LCD screen 5", with adjustable brightness,
  - Keypad lock,
  - Air temperature control mode temperature range: 20 to 37°C,
  - Temperature range of skin temperature control mode: 34 to 37°C,
  - Ability to monitor two skin temperatures simultaneously,
  - Noise level under the dome < 47 dBa (without oxygen servo control),
  - Selectable menu in PL, EN, FR, IT, DE, ES, NN, HU, DA, RU, TR, FI, PT, JA, ZH,
- Included:
  - User manual in Polish (PDF),
  - DRAGER body temperature sensor,
  - Condensation jar with a capacity of 800 ml,
  - Power cord,
- Dimensions:
  - Mattress: 36.5 x 73 x 2.3 cm,
  - Chamber: 40 x 83 x 50 cm,
  - Overall: 75 x 104 x 132.6 / 152.7 cm,
- Weight: 93 kg,
- Cat 02,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN DUTCH  
Мы говорим по русски  
Мы говорим по украински

