



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER ISOLETTE C2000 Incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 300 W,
- DRAGER Isolette C2000 is a neonatal incubator that provides a high level of thermoregulation and an ideal micro-environment for the development of newborns,
- Equipped with a new medical approved mattress in cream color,
- Has built-in systems:
 - Maintaining constant temperature,
 - Air humidity regulation,
 - Alarm,
- Equipped with:
 - 6 nursing openings for easy access to the patient,
 - Distilled water container,
 - Pull-out tray,
 - 1 storage compartment,
 - Rails for attaching additional accessories,
 - Grommet seal for tubes and drains - 10 pieces,
 - Socket for connecting two skin temperature sensors (no sensors included),
 - Foiled control panel,
 - Selectable menu in the following languages: PL, NR, PA, PT, JP, IT, SP, DE, FR, EN, SW, FN, GR, CZ , SL,
 - Power sockets for connecting additional devices - 6 pcs,
 - Built-in scale,
- Specifications:
 - Infinitely variable manual change of the recliner angle to Trendelenburg and anti-Trendelenburg positions within 12º,
 - Possibility to adjust the humidity of the air in the chamber in the range: 30 - 95 %,
 - Air temperature: 20 - 37º C,
 - Skin temperature: 34 - 37 º C,
 - Pressure: 40 - 150 psi,
 - Oxygen: 20 - 65%,
 - Heating time: < 35 min,
 - Fresh air flow rate: 30 L / min,
 - Noise level inside the chamber: < 47 dB,
 - Mobile on 4 wheels with Ø 12 cm, all with brakes,
- Included:
 - User manual in English (PDF),
 - Drager body temperature sensor,
 - Power cord,
- Dimensions:
 - Mattress: 37 x 72 x 2 cm,
 - Chamber: 45.5 x 83.5 x 49 cm,
 - Overall: 70 x 105 x 140 cm,
- Weight: 79kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

