



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER ISOLETTE 8000 Neonatal incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in USA, (2012),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Isolette 8000 enables continuous monitoring of central and peripheral body temperature of newborns, is an ideal incubator for neonatal intensive care units,
- Medical approved mattress in white color,
- Advantages:
 - Dual air curtains and advanced servo-regulation system maintain stable conditions inside the incubator during patient access,
 - The thermo-monitoring system gives you the opportunity to observe the baby's central and peripheral skin temperature, also in servo control mode,
 - Trends are clearly presented on the LCD display.
 - Easy access to the newborn,
 - Convenient electronic height adjustment of the incubator,
 - Manual adjustment of the chamber angle with a knob (Trendelenburg or anti-Trendelenburg position) within the range: $\pm 12^\circ$,
 - Low weight, compact design,
- Specifications:
 - LCD screen 5", with adjustable brightness,
 - Air temperature control mode temperature range: 20 to 37 $^\circ$,
 - Temperature range of skin temperature control mode: 34 to 37 $^\circ$ C,
 - Ability to monitor two skin temperatures simultaneously,
- Equipped with:
 - Swiveling front and rear access panel,
 - 6 nursing holes for easy access to the patient,
 - 10 x grommet seals for wires and drains,
 - Socket for connecting two skin temperature sensors,
 - Rails for attaching additional accessories,
 - 3 storage drawers,
 - Self-adjusting casters, and ergonomic foot brakes allow the incubator to be moved easily, even in limited ward space,
- Noise level under the dome < 47 dBA (without oxygen servo control),
- Selectable menus in EN, PL, DE, IT, ES, JA, DE, DK, NORSE, POL, PT, SWE, FI, GRE, CS, SK, RU, HU, TR,
- Included:
 - User manual in Polish (PDF),
 - Power cord,
 - Condensate jar with a capacity of 800 ml,
 - 1 x patient temperature sensor,
 - Hose for oxygen administration;
- Dimensions:
 - Mattress: 36.5 x 73 x 2.3 cm,
 - Chamber: 40 x 83 x 50 cm,
 - Overall: 75 x 104 x 132.6 / 152.7 cm,
- Weight:
 - 93 kg,
- Cat 03,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

