



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KNÜRR Endoscopy cart (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Specifications:
 - Top drawer lockable,
 - Power switch on the front of the cart,
 - Equipped with electrical installation,
 - Endoscope holder with height adjustment in range: 60 cm,
 - Removable rear housing,
 - Mobile on 4 casters including 2 with brake,
- Included:
 - Power cords (4 pcs),
- Dimensions:
 - Top: 42 x 44 cm,
 - Shelves: 45 x 45 cm,
 - Top drawer: 34 x 37cm,
 - Lower drawer: 36 x 38cm,
 - Overall dimensions: 60 x 55 x 122 cm,
- Weight: 49kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

