



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER ISOLETTE C2000 Newborn incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 300 W,
- DRAGER Isolette C2000 is a neonatal incubator that provides a high level of thermoregulation and an ideal micro-environment for the development of newborns,
- Equipped with a new medical approved mattress in cream color,
- Has built-in systems:
 - Maintaining constant temperature,
 - Air humidity regulation,
 - Alarm,
- Equipped with:
 - 6 nursing openings for easy access to the patient,
 - Pull-out tray,
 - Distilled water container,
 - Rails for attaching additional accessories,
 - Grommet seal for tubes and drains - 10 pieces,
 - Socket for connecting two skin temperature sensors (no sensor included - can be retrofitted at additional cost),
 - Power sockets for connecting additional devices - 6 pcs,
 - Built-in patient weight (accuracy in 10 gram increments) - measurement up to a maximum of 7 kg,
 - Foiled control panel, has a small scratch on the screen, however, it does not significantly affect the comfort of operation of the device (photo),
- Specifications:
 - Infinite manual change of the recliner angle to Trendelenburg and anti-Trendelenburg positions within 12°,
 - Air temperature: 20 - 37° C,
 - Skin temperature: 34 - 37 ° C,
 - Pressure: 40 - 150 psi,
 - Heating time: < 35 min,
 - Fresh air flow rate: 30 L / min,
 - Oxygen: 20 - 65%,
 - Noise level inside the chamber: < 47 dB,
 - Selectable menu in language: PL, NR, PA, PT, JP, IT, SP, DE, FR, EN, SW, FN, GR, CZ , SL,
- Mobile on 4 wheels with Ø 12 cm all with brakes,
- Includes:
 - User manual in English (PDF),
 - Power cord,
- Dimensions:
 - Mattress: 37 x 72 x 2 cm,
 - Chamber: 45.5 x 83.5 x 49 cm,
 - Overall: 70 x 105 x 140 cm,
- Weight: 79 kg,
- Cat 03,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

