



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHYSIO-CONTROL LIFEPAK 15 Defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good ,
- Visual condition: very good,
- Real photos of the product,
- Made in USA, (2014),
- Power supply: 230 V,
- Through 2 built-in batteries: 10.8 V / 6 Ah / 65 Wh - can be exchanged for new ones at additional cost,
- Two-phase,
- Portable, equipped with:
 - Cushioning handle for carrying,
 - Screen with a double layer of protection to withstand heavy impacts and drops,
 - LCD color liquid crystal display with diagonal: 8,4",
 - English language control panel,
 - Defibrillator menu in English,
 - Built-in printer with paper width: 10.5 cm,
- Wireless connection with other devices via Bluetooth,
- Energy range in manual mode: 2 - 360 joules,
- Built-in memory capacity:
 - 360 minutes of continuous ECG monitoring,
 - 90 minutes of continuous data monitoring on all channels,
 - 400 individual parameter curve events,
- Six modes of operation:
 - AED mode - for automatic ECG analysis and guided protocol for patient and cardiac arrest management,
 - Manual Mode - for manual defibrillation, synchronized cardioversion, non-invasive electrostimulation, and vital signs and ECG monitoring,
 - Archive mode - to access stored patient data,
 - Setup Mode - for changing default settings of operational functions,
 - Service Mode - intended for authorized personnel and used to perform diagnostic tests and calibrations,
 - Demo mode - for simulating curves and change graphs, for demonstration purposes,
- Features:
 - SYNC - Synchronization - Cardioversion,
 - PACER - Stimulator (activates the pacing function),
 - ANALYZE - activates the SAS (Shock Advisory System),
 - TRANSMIT - transmits patient data,
- Equipped with the following modules and cabling:
 - ECG - Main cable for 10-lead ECG,
 - SpO2 - Cable with finger clip sensor,
 - NIBP - Reusable cuff,
 - CO2 - Sensor including filter,
 - QUIK - COMBO electrode adapter with electrode included,
- Selectable menu in language: EN, DE,
- Keypad on the control panel in German,
- Included:
 - User manual in language: EN, EN (PDF),
 - Hard spoons,
 - Bags for defibrillator accessories,
 - Batteries: 2 x 10.8 V / 6 Ah / 65 Wh,
 - Power supply,
 - Tester,
 - Wiring,
 - Power cord,
- Dimensions of the defibrillator: 23 x 40 x 34 cm,
- Dimensions with bag: 22 x 64 x 34 cm,
- Weight of the set: 13 kg,
- Cat. D126,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

