



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CRYO - T ELEPHANT MINI (50l) Apparatus (New)**

- New,
- Made in Poland,
- CRYO-T Elephant Mini - A proven, economical solution for topical cryotherapy for SPA centers and rehabilitation clinics,
- Quick start - Quick START allows you to start working right after turning on the device. In a short period of time, the device generates the cold between -110°C -160°C required for effective treatment. The desired temperature is reached 1 minute after pressing the START button,
- Patented technology - High-tech cryogenic tubing (protected by international patent) provides full freedom of movement. It does not stiffen during the procedure, as a result, the operator's hands do not get tired. The insulating material used in the cryogenic cable prevents it from freezing,
- Savings of up to 40% - The device is available with specially developed heater technology, which provides high savings in refrigerant consumption of up to 40%,
- Non-freezing hose - patented technology - The non-freezing, patented cryogenic hose guarantees unprecedented comfort of use and full freedom of movement,
- Up to -160° C - Thanks to the use of nitrogen as a cooling agent, the treatment temperature is as low as -160° C,
- Fast and safe treatment - Treatment time in one area of the body ranges from 30 seconds to 5 minutes, on the face from 2 to 3 minutes. Patients with increased fat or muscle tissue may require extended treatment up to 10 minutes. The proper application of cold lies in the experience and knowledge of the doctor/therapist, as well as careful and constant observation of the reaction of the tissues subjected to cold,
- Application:
  - Physiotherapy and Rehabilitation,
  - Aesthetic Medicine,
  - Rehabilitation,
  - Sports Medicine,
  - Wellness,
- Technical data:
  - Power consumption: 230V,
  - Refrigerant: LN2,
  - Cold jet temperature: up to -160°C,
  - Adjustment of jet intensity: 5 levels,
  - Available tank capacity: 50 l, 35 l, 11 l,
  - Coolant level indicator: built-in,
  - Nozzle: traditional (without indicator and skin temperature measurement),
- Included:
  - Apparatus,
  - 35l tank,
- Has a current Technical Passport issued,
- Warranty 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

