



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



ANETIC AID QA3 Transport stretcher (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in UK,
- New upholstery with medical approval in black,
- Adjustable in many planes (hydraulically and through gas springs):
 - Bed height in the range: 67 - 92 cm
 - Trendelenburg position in the range: 14°,
 - Backrest angle using in the range of 0° to 90°,
- Features:
 - X-RAY translucency,
 - Folding handrails,
 - Retractable IV stand in the handrail, Adjustable from: 140 - 180 cm,
 - Plastic base,
 - Round rubber bumpers on the corners,
 - Mobile on four casters with diameter: 15 cm,
 - Central brake and the possibility of directional travel by means of the extendable fifth wheel,
- Dimensions:
 - Overall: 77 x 209 x 84.5/114 cm,
 - Mattress: 54 x 190 x 8 cm,
- Maximum load capacity: 160 kg,
- Weight: 125 kg,
- Cat. 02,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

