



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **KART-MAP Office cabinet (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Office closet with doors that swing outward and lock. Built-in storage compartment in the upper part of the cabinet with doors opening and locking with a key. Shelves can be adjusted in height. Solid workmanship and stable construction,
- Includes:
  - 3x shelf with fasteners,
  - 2x Key for locks,
- Dimensions:
  - Total: 90 x 42 x 192 cm,
  - Shelf: 90 x 35 cm,
  - Storage compartment: 90 x 35 cm,
- Weight: 60kg,
- Cat 01,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

