



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**SONOSCAPE S2 Ultrasound scanner (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in PRC, (2017),
- Power supply: 230 V,
- Emergency power supply: via built-in 14.8V, 10.4 Ah battery (can be replaced with a new one at additional cost),
- Frequency: 50/60 Hz,
- Power: 10 A,
- SONOSCAPE S2 is a state-of-the-art ultrasound machine with amazing image quality, built on the basis of the latest solutions. This unit uses the new microScan multibeam imaging technique, which dramatically improves image quality and clarity. The S2 has a powerful Doppler application, enabling comprehensive "in color" diagnostics. Thanks to the extraordinary sensitivity of the color Doppler, examining even small flows becomes very easy. It features unparalleled functionality, a large 15-inch LCD monitor, programmable keys, intuitive menus and an examination archiving system with a patient database. The SonoScape S2 works with a wide range of broadband heads for various types of examinations.
- Specifications:
  - 15" monitor,
  - Micro Convex C611 head,
  - Two active probe ports,
  - Adjustable probe cable holder,
  - Sockets:
    - ECG,
    - VIDEO OUT,
    - FOOT SW,
    - PRINT,
    - VGA-OUT,
    - ETHERNET,
    - USB x2,
    - S-VIDEO,
  - Mobile, functional table equipped with:
    - Operator handle
    - 4x head mount,
    - 2x shelf,
    - Accessory storage compartment,
    - Adjustable head cable holder,
    - 4x braking casters,
- Features:
  - Spatial Compound Imaging - simultaneous scanning of the ultrasound beam from multiple angles (compound type) so-called crossed ultrasound imaging,
  - Tissue Harmonic Imaging (THI) - broadband harmonic imaging,
  - Color Doppler - new generation dynamic wide-band Color Doppler mode provides accurate analysis of blood flow morphology with excellent sensitivity and exceptional quality,
  - Panoramic Imaging- optional panoramic imaging,
  - Real-time 4D Imaging- real-time 3D imaging using a volumetric probe (optional),
  - Examination modes: abdominal, obstetrics, gynecology, small organ, musculoskeletal, rheumatology, vascular, urology, cardiology, etc,
  - Clip-board function - quick and easy review of images during examination,
  - M-tuning function - Intelligent image optimization at the touch of a button,
  - Selectable menu in PL, EN, CN, RU, DE, IT, FR, ES, RU, HU, PA,
- Included:
  - User manual in: Polish, English (PDF),
  - Power cord,
  - Micro Convex C611 head,
  - Carrying case with handle,
  - Ground cable,
  - S-VIDEO cable,
  - Mobile cart,
- Dimensions:
  - USG: 35 x 15 x 42 cm,
  - Cart: 52 x 62 x 90 cm,
  - Cart top shelf: 27 x 14 cm,
  - Cart bottom shelf: 38 x 22 cm,
- Total weight: 41.5 kg,
- Cart: 25 kg,
- USG: 9.5 kg,
- Accessories + instrumentation: 7 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

