



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC BIPHASIC LIFEPAK 12 Defibrillator (Cat.78) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- **Real product photos,**
- Made in the USA,
- Power supply: two replaceable 12V 1.2 Ah batteries (**one included, the second one can be added for an additional fee**)
- Frequency: 50/60 Hz,
- Power: 180 W,
- **Possibility to select menu language: EN, FR, DE, ES, SW, IT, CN, NL, PL, DK, NR,**
- Keyboard on the control panel in: English,
- Built-in printer for paper width: 11 cm,
- Equipped with a carrying handle,
- **Features:**
 - SYNC - Synchronization - Cardioversion,
 - ADVISORY - INFO (advice) - Activates the continuous patient surveillance system (CPSS),
 - ANALYZE - Activates the SAS (Shock Advisory System),
- **Equipped with the following modules with wiring and sensors:**
 - ECG - 10-lead cable,
 - SPO2,
 - Cuff with hose,
- **For an additional fee, the defibrillator can be equipped with:**
 - Hard spoons,
 - Additional battery,
 - External charger from AD Elektronik visible in the last photo,
- **Included in the set:**
 - User manual in Polish (PDF),
 - Quick-Combo therapy cable,
 - Removable screen and operator keyboard cover,
 - Retrofitted transport bags for equipment,
- Dimensions:
 - Defibrillator: 20 x 40 x 34.5 cm,
 - With bag: 23 x 62 x 34.5 cm,
- Weight: 8.5 kg,
- CAT. 78,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,
-

