



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



PHILIPS INTELLIVUE MP50 Vital function and anesthetic gas monitor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- IntelliVue MP50 is a cardiac monitor that monitors multiple vital functions,
- It has:
 - Panel screen, made of TFT LCD technology with a diagonal of 12" and working in SVGA resolution, has minor scratches that do not affect the operation of the patient monitor,
 - Operation by touch screen, and navigation joystick allows direct access to many functions with simple commands. Up to 6 curves can be displayed on the screen, while 12 ECG leads can be displayed on the screen,
 - The convenient handle and sturdy case make it easy to carry,
 - Custom display options allow you to view and analyze data in graphical and numerical formats, compile real-time measurements and trends, and freely configure each screen element - from curves to parameter names,
 - IntelliVue technology - access information from various systems in the hospital network - pharmacy, laboratory, ECG management, radiology department PACS, server-based hospital systems and others - and deliver it to the patient monitor,
 - A navigation dial makes it easy to enter information and navigate the screen,
 - Audible and visual alarms,
 - Defib Synch socket, LAN,
- Equipped with PHILIPS Model M3001A module with wiring:
 - ECG /RESP - Brand new 3 - lead cable with extension,
 - SPO2 - New silicone sensor,
 - NBP - Cuff with air tubing,
 - TEMP - Temperature connector socket, (no sensor included, can be retrofitted at additional cost),
- PHILIPS IntelliVue gas analysis module, model G5 M-1019A:
 - Continuous real-time monitoring of respiratory gases and anesthetics (simultaneously),
 - IntelliVue G5 integrates completely with Philips IntelliVue MP40/50/60/70/80/90 anesthetic patient monitors,
 - Automatic identification of up to two anesthetic agents,
 - Inspiratory and expiratory levels of all gases,
 - Respiratory rate (awRR) calculated from the CO2 curve,
 - Minimum alveolar concentrations (MAC),
 - Long-term stability enables simple, infrequent calibration,
 - Individual alarm thresholds for each parameter,
 - Quick access to the most commonly used functions,
- English language menu,
- Included:
 - User manual in Polish (PDF),
 - Power cord,
- Dimensions: 37 x 20 x 43 cm,
- Weight: 11 kg,
- Cat 04,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

