



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DATEX-OHMEDA ISOTEC 5 Isoflurane vaporizer (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in England,
- Power supply: 230 V,
- GE Datex-Ohmeda TEC 5 Isoflurane is designed for off-circuit use in continuous flow techniques for inhalational anesthesia. The TEC 5 is temperature-, flow- and pressure-compensated, so its performance remains relatively constant despite evaporation-induced cooling,
- The TEC 5 evaporator can be installed in Datex-Ohmeda anesthesia machines equipped with the Selectatec attachment manifold. No tools are needed to perform this operation, and the anesthesia machine does not need to be taken out of service,
- Ability to read the liquid anesthetic level thanks to an indicator on the front of the evaporator,
- Liquid capacity (dry wick): 300 ml,
- Liquid retained by the wick: 75 ml,
- Temperature operating range: 18°C to 35°C,
- Dimensions: 6.5 x 11 x 25 cm,
- Weight: 7 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

