



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NANASONIC Ophthalmic glass drilling machine (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 15 W,
- Opposite taper drill system for simultaneous drilling from both sides for precise cutting,
- Fine adjustment for precise drilling,
- Technical specifications:
 - 1250 RPM,
 - Height adjustable slide stand,
 - Side-to-side adjustable slide stand,
 - Manual lever for adjusting the extension of the needle,
- Included:
 - Power cord,
 - Gasket,
 - 2x Allen,
 - 1x slide support,
 - 1x container with funnel for cooling substance,
 - Adjustment extension,
- Dimensions: 20 x 20 x 30 cm,
- Weight: 4 kg,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

