



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**MALEVESTIO VIVO Hospital bed (ICU) (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Italy, (2018),
- Power supply: 230 V,
- Emergency power supply: through built-in batteries,
- Frequency: 50/60 Hz,
- Power: 750 VA,
- VIVO is an electric bed for intensive care, with 4 movable segments, with height adjustment by means of high-performance telescopic columns. The bed's recliner is equipped with a patented, INNOVATIVE PAS SYSTEM (to prevent the patient from slipping), which COMPENSATES the patient's pushing towards the legs of the bed. The leg-side top has an integrated multifunctional touchscreen to control bed adjustment and the VI-MAT air mattress pumping system and VI-WEIGHT patient weighing system,
- The air mattress is pumped using a compressor built into the leg-side of the bed, making it easy to clean and reducing the overall size of the air supply system for the mattress. With the side-tilt function, rotational treatment can be programmed for patients with pulmonary complications. All controls are ergonomically arranged and labeled with intuitive icons easily understood by both medical personnel and the patient,
- In addition to the digital touchscreen, the bed can be operated via control panels located on the outside of the side rails and equipped with LED lighting for better visibility at night. The LED flashes in the event of an alarm or nurse call. Patient button functions, located on the inside of the side rails, can be interlocked by nursing staff.
- Side tilting can be controlled by pedals, making it easier for the operator to care for the patient. The control panel also allows programming of rotational therapy functions.
- Technical specifications:
  - External operator panels on both sides that control bed functions,
  - Integrated patient panel with intuitive icons,
  - Touch operator panel on the patient's foot side to manage all functions and bed settings,
  - Electric backrest adjustment: 0 / + 68°,
  - Trendelenburg position 16°,
  - Anti-Trendelenburg position 15°,
  - Electric tilting of bilateral side tilt: 15 / - 15°,
  - Electric adjustment of knee break: 0 / + 25° (incrementally, relative to the angle of tilt of the central segment),
  - Thigh segment angle: 0 / - 55° (relative to the angle of the knee rests) ,
  - Ability to electrically extend the bed within 23 cm,
  - Height adjustment within the range: 48 - 82 cm,
  - Touch-reset button,
  - Bilateral CPR lever,
  - Bilateral CPR pedal,
  - X-ray translucency,
  - Integrated patient weight,
  - Foot control,
  - Built-in mattress,
  - Anti-slip trowel system
  - Patient anti-slip system,
  - Fifth wheel,
  - Central braking system with built-in wheel alarm,
- Includes:
  - Power cord,
  - Spare mattress cover,
  - Low-pressure air mattress, with dual filling mode.
  - Replacement cover,
- Dimensions:
- Overall: 108 x 219 cm,
- Platform length: 198 cm,
- Platform width: 90 cm,
- Clearance under the bed: 13 cm,
- Maximum load capacity:
  - Patient: 190 kg,
  - With working load: 260 kg,
- Total weight: 240 kg,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

