



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE IES 300 Smoke evacuator (Reconditioned)

- Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Real photos of the product,
 - Made in Germany,
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - Power: 3.0 A
 - The ERBE IES 300 smoke evacuator extracts harmful volatile fractions from the surgical field and increases occupational safety by minimizing the inhalation of unwanted substances by medical personnel,
 - Programmed and controlled by the diathermy, but manual control of the unit's operation via a control panel is also possible,
 - Low noise level during operation despite high flow rate,
 - Holistic safety concept: e.g. automatic interlock in case of tissue aspiration (at high flow rate),
 - Long filter life (35 h) at highest flow rate,
 - Air flow rate (free flow): 1600 l/min,
 - Maximum flow rate: 650 l/min,
 - No main filter (can be retrofitted at additional cost),
 - Includes:
 - User manual in English (PDF),
 - Power cord,
 - Diathermy-evacuator connector,
 - Dimensions: 41 x 37 x 21 cm,
 - Weight: 14 kg,
 - Has an up-to-date inspection and is ready for use,
 - Issued Technical Passport (Service Report) valid for 12 months,
 - Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
 - Financing options (Poland only): Installments, Leasing, Loan,
- If you have any questions, please do not hesitate to contact us!
- In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

