



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **FERNO COMPACT Cardiac chair (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in England,
- Safety belts with Easy-Clip buckle and raised seat panel increase patient comfort and safety,
- The chair is foldable, compact and lightweight,
- Easy to clean,
- Seat and backrest made of nylon, flame retardant material,
- Stainless steel construction,
- Has two wheels and two lifting loops,
- Wheels set inside the frame to a minimum,
- Dimensions:
  - Height (folded): 93.5 (77.5)cm
  - Width (folded): 44.5 (44.5)cm
  - Depth (folded): 42 (28.5)cm
- Maximum load capacity: 200 kg,
- Weight: 9.1 kg,
- It has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

