



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **ESPE ELIPAR CLASSIC Polymerization lamp (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- New 12.5 V 75 W bulb,
- Polymerization wavelength range: 400 -500nm,
- Polymerization in two-stage mode,
- Technical specifications:
  - Polymerization wavelength range: 400 -500nm,
  - Polymerization in two-stage mode,
  - Light intensity: 800 nW/cm<sup>2</sup> (average value),
- Included:
  - User manual in Polish (PDF),
  - Power cord,
- Total dimensions of the fixture:
  - 43 x 28 x 15 cm,
- Weight: 2.5 kg,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



