



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FRASENIUS KABİ APPLIX SMART D Enteral feeding pump (Reconditioned)

- Used, very good technical condition,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: Through new built-in 4.8 V batteries,
- Frequency: 50/60 Hz,
- Power: 10 VA,
- The pump has a 45 x 60 mm display,
- Keypad on the control panel in English,
- The Applix Smart D pump is designed exclusively for enteral feeding,
- It is very easy to use and can be used both by home patients and in hospital treatment,
- The pump has a continuous feeding program and functions to ensure patient safety,
- Accuracy of feeding product delivery at a rate of more than 3 ml/h,
- Speed of feeding product delivery: 1 - 600 ml / h;
- Volume of nutritional product to be administered: 1 - 5000 ml;
- Feed delivery accuracy: 10 %;
- Maximum operating pressure: 200 kPa;
- Battery-powered pump operating time at 125ml / h: 24 h;
- Easy cleaning with removable door,
- Adjustable alarm volume,
- The pump can be mounted universally, vertically and horizontally to round and square rails,
- Equipment:
 - Power cord,
 - User manual in Polish (PDF),
 - Pump holder,
 - Docking station,
- Can be retrofitted with a soft-pack enteral diet kit at an additional cost,
- Dimensions: 21 x 15.5 x 14 cm,
- Weight: 1.1 kg,
- Has a current Technical Passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 1 month for the international market,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

