



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TOPCON RM-A6000 Autorefractometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 130 VA,
- Topcon RM - A6000 is an autorefractometer for computerized examination and evaluation of visual acuity by taking measurements of the curvature of the cornea and lens of the eye,
- The autorefractometer test itself consists of injecting a beam of infrared radiation into the eye and analyzing the reflection of the beam off the structures of the eye. The computer built into the device performs a thorough analysis of the reflections obtained, based on which it generates data showing the patient's specific refractive defect,
- Product features:
 - Manual height adjustment:
 - Podium via joystick,
 - Chin support by means of a screw,
 - Practical joystick makes it easy to set parameters for the test,
 - Momentary button starts the selected test,
 - 4" TV display,
 - Adjustable brightness and contrast,
 - Built-in printer for 6 cm wide thermal paper,
 - RS - 232C connector for data transmission to computer,
 - Possibility of displaying 6 -five optotypes for eye examination,
 - Switching between optotypes using foil buttons above the TV screen,
- Functions available from the control panel:
 - Switching between manual / automatic mode,
 - Zoom in / out button,
 - Print button,
- Parameters accessible from the hidden panel under the flap:
 - Patient VD: 12 or 13.75 mm,
 - Measurement step (STEP): 0.12 or 0.25,
 - Cylinder (CYL): +, - or MIX,
- English language menu,
- Included:
 - Power cord,
- Overall dimensions: 27.5 x 45 x 47.5 cm,
- Weight: 21 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



PRAISTON



PRAISTON



PRAISTON



PRAISTON



PRAISTON



PRAISTON



PRAISTON