



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment and cosmetic chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- Real photos of the product,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Upholstery with medical certificate in white color (it is possible to replace the upholstery with a new one from our color pattern for an additional fee),
- Ideal for classic cosmetic procedures,
- Electric up/down adjustment with remote control in the range: 64 - 84 cm,
- Armrests automatically adjustable along with the angle of the backrest,
- Manual adjustment in the range:
 - Backrest: 5° / 55°,
 - Footrest: 0° / - 34°,
 - Removable headrest,
- Included:
 - Power cord,
 - Up/down control remote control,
- Dimensions:
 - Headrest: 40 x 25 x 6 cm,
 - Backrest: 61 x 70 x 11 cm,
 - Seat: 61 x 51 x 11 cm,
 - Footrest: 60 x 62 x 12 cm,
 - Armrests: 10 x 73 x 6 cm,
 - Length in reclining position: 185 cm,
 - Width: 85 cm,
 - Height: 125/137 cm,
 - Base: 57 x 106.5 x 24 cm,
- Weight: 65 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

