



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS UCR Insufflator (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 40 VA,
- OLYMPUS UCR insufflator increases patient comfort during and after colonoscopy,
- Carbon dioxide is absorbed by the body 150 times faster than air,
- Easy one-button operation,
- Pressure level display,
- Three controlled airflow settings available,
- Max, supply pressure: 45 kPa,
- Can be retrofitted at an additional cost with the visible in the last photo:
 - Water bottle (Olympus MAJ-902),
 - Gas tube (Olympus MAJ-1741),
- Included:
 - CO2 bottle hose,
 - Power cord,
- Dimensions: 12.5 x 30 x 15 cm,
- Weight: 5 kg,
- Comes with a valid Technical Passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

