



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MINDRAY VS-800 Patient monitor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in PRC,
- Power supply: 230 V,
- Emergency power supply: through 1 new built-in 12 V 2.3 Ah battery,
- Frequency: 50/60 Hz,
- Mindray VS-800 is a multi-parameter patient monitor for adult, pediatric and neonatal patients. Parameters that can be monitored include: SpO2, PR, NIBP and TEMP. The VS -800 has a large LCD screen that can display high-resolution waveforms. The monitor can store up to 100 patient IDs and 1,200 records, and can transfer data to a PC or central station for review. The Mindray VS-800 is equipped with a removable and rechargeable lithium-ion battery that provides up to 10 hours of continuous operation.
- Features:
 - Suitable for adult, pediatric and neonatal patients,
 - Clinically proven SpO2 technology,
 - SmarTemp thermometer technology,
 - Adjustable audible and visual alarms,
 - Central station network,
 - Removable and rechargeable battery,
- Technical specifications:
 - Display:
 - 7-segment LEDs: systolic pressure, diastolic pressure, mean pressure, SpO2, temperature and heart rate,
 - 3.4-inch FSTN LCD display: resolution: 320 x 160, 0.24 x 0.24 mm spacing,
 - Trace: 1 plethysmogram curve,
 - Sweep speed: 25 mm/sec,
 - NIBP:
 - Method: Automatic oscillometric,
 - Range: Adults: 10-270 mmHg, Children: 10-200 mmHg, Newborns: 10-135 mmHg,
 - Overpressure protection: yes,
 - Heart rate: 40-240 beats per minute,
 - Heart rate resolution: 1 bpm,
 - Heart rate accuracy: ±2 beats per minute,
 - SPO2:
 - Range: 0-100%,
 - Resolution: 1%,
 - Accuracy:
 - Mindray SpO2: adults/children: ±2% (70-100%), neonates ±3% (70-100%),
 - Masimo SpO2: adults/children: ±2% (70-100%), neonates ±3% (70-100%),
 - Nellcor SpO2: accuracy: ±2-3.5% (70-100%),
 - Heart rate:
 - Mindray SpO2: range: 20-254 beats per minute,
 - Masimo SpO2: range: 20-240 beats per minute,
 - Nellcor SpO2: range: 20-250 beats per minute,
 - Resolution: 1 bpm,
 - Accuracy: ±3 beats per minute,
 - Temperature sensor:
 - Range: 25-44°C (77-112.2°F),
 - Resolution: In monitoring mode: 0.2°F (0.1°C),
 - Accuracy: In monitoring mode: 77-89.6°F (25-32°C): ±0.3°F (±0.2°C); 32-44°C: ±0.2°F (±0.1°C),
 - Measurement time: <12s,
 - Recorder:
 - Type: Built-in, Thermal chip,
 - Record width: 48 mm (50 mm registration paper),
 - Print speed: 25 mm/sec,
- Keypad on control panel in English,
- Included:
 - Power cord,
 - NIBP cuff + cable,
 - MINDRAY SPO2 sensor,
- Dimensions: 17 x 17 x 24 cm,
- Weight: 3.5 kg
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

