



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **WOSON HSV IID7 Autoclave (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in PRC. (2011 r.),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Rated power: 1400 VA,
- Device class: B,
- The unit is designed to sterilize all packaged and unpackaged, solid, concave type A and porous items,
- Clear digital LCD display with automatic shutdown to save energy,
- Operating programs:
  - Rapid sterilization program (PROMPT) - temperature 134°C; pressure 2.1 bar; sterilization program time 3 min; drying time 10 min,
  - Instrument set sterilization program (PACKAGE) - temperature 134°C; pressure 2.1 bar; sterilization program time 5 min; drying time 10 min,
  - Sterilization program for cotton or selected instruments (COTTON) - temperature 134°C; pressure 2.1 bar; sterilization program time 15 min; drying time 17 min,
  - Sterilization program for rubber or selected instruments (RUBBER) - temperature 121°C; pressure 2.1 bar; sterilization program time 15 min; drying time 17 min,
- Tests:
  - Vacuum Test,
  - Bowie&Dick Test,
  - Helix Test,
- Technical specifications:
  - Autoclave chamber made of stainless steel,
  - Chamber capacity: 7 liters,
  - Maximum operating pressure: 2.5 bar,
  - Minimum working pressure: - 0.9 bar,
  - Maximum temperature: 145 °C,
  - Volume: < 50 dB,
  - Built-in input for printer connection,
- Includes:
  - User manual in Polish (PDF),
  - Power cord,
  - Sterilization trays - 3 pieces,
  - Tray holder,
  - Water supply tube,
- Dimensions:
  - Trays: 13 x 26 x 2 cm,
  - Chamber: Ø 17 x 32 cm,
  - Overall: : 37.5 x 55 x 38 cm,
- Weight of set: 45.5 kg,
- It has a current inspection and is ready for operation,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

