



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **KTG BIONICS IFM-500 KTG unit (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Korea,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- IFM-500 is a mobile and lightweight monitor, used in obstetrics and gynecology for prenatal diagnosis of the fetus,
- It is designed to monitor and record the heart rate of 2 fetuses and uterine contractions,
- The patient-operated event marker and the ability to monitor fetal movements make the examination painless, and the alarm system immediately alerts the operator to any changes,
- Equipped with:
  - Digital displays showing:
    - Accurate and stable FHR (fetal HeartRate),
    - UC (uterine contraction),
    - FM (fetal movement) extracted automatically,
  - Built-in thermal printer that automatically prints data for periods of 10, 20, 30, 40, 50, 60 minutes, so the operator can get complete information about the patient's condition,
  - Convenient carrying handle,
- Specifications:
  - Operating frequency: 2.0 MHz,
  - Intensity: 10 mW/cm<sup>2</sup> or less,
  - FHR measurement:
    - Dual display and FHR recording,
    - Measurement range: 50-240 bpm (beats per minute),
    - Fetal heartbeat sound,
    - Fetal heart rate light,
    - Signal quality indicator,
    - FHR alarm function,
    - Automatic recognition of the transducer used,
  - FM measurement:
    - Automatic measurement by Doppler shift frequency,
  - UC measurement:
    - External measurement (guard ring method),
    - Calibration function,
    - Display of UC values on LED strips,
    - Displaying and recording UC values,
  - Printer:
    - Thermal 8 dot/mm,
    - Print speed: 1.2 or 3 cm/min,
    - Print contrast: 1 (normal), 2, 3, 4,
    - Auto. prints 10, 20, 30, 40, 50, 60 min,
- Control panel keyboard in English,
- Included:
  - User manual in English (PDF),
  - Power cord,
  - DOP head - 2 pcs,
  - UC head - 1 pc,
  - Patient button,
  - Head attachment straps - 3 pcs,
- Dimensions:
  - Monitor: 22 x 23 x 22 cm,
- Weight: 12.5 kg,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

