



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCAN MODULE SYSTEM Medical trolley (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the Netherlands,
- Power supply: 230 V,
- SCAN MODULE SYSTEM medical trolley on 4 casters (2 braking),
- Pull-out drawer with key locking option,
- Built-in handles with the ability to move the assistant comfortably,
- Construction made of aluminum and plastic,
- Included:
 - Drawer keys (4 pcs),
- Dimensions:
 - Top shelf: 46 x 70 cm,
 - Middle shelf: 40 x 60 cm,
 - Lower shelf: 40 x 60 cm,
 - Pull-out drawer: 34.5 x 54 x 20 cm
 - Total: 46 x 70 x 90 cm,
- Weight: 45 kg,
- Has an up-to-date inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

