



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **MPW MPW-350R Laboratory centrifuge (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 870 W,
- MPW-350R centrifuge is used in medical, biochemical laboratories and food industry,
- Features of MPW-350R:
  - Microprocessor control,
  - Choice of 99 memory programs,
  - Selection of 10 acceleration and deceleration curves,
  - Automatic rotor identification,
  - Fault identification system,
  - Registration of operating parameters - RS 232 (optional),
  - Eco-friendly and efficient cooling system,
  - Pre-cooling capability,
  - Product compliance with international safety conditions,
- Technical specifications:
  - RPM range: 300 - 15000 / rpm,
  - Spinning time: 0 - 99 min,
  - Temperature control range: -10 / + 40 °C,
- Included:
  - User manual in Polish (PDF),
  - Power cord,
- No rotor included - it is possible to retrofit the centrifuge with a rotor and selected accessories (cups/adapters/samples) for an additional fee,
- Dimensions: 41 x 64 x 37 cm,
- Weight: 59 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

