



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



UNIPHY PHYACTION PERFORMA Short-wave diathermy (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the Netherlands,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 64 W,
- UNIPHY PHYACTION PERFORMA is a short-wave diathermy used to raise tissue temperature and exert therapeutic effect,
- Easy and intuitive operation,
- Applicators placed on stable, articulated masts adjustable in multiple planes,
- Even heating makes treatment very comfortable for patients,
- Technical specifications:
 - 2 input channels,
 - 2 applicators:
 - 100 W,
 - 200 W,
 - Emitting wave frequency (F): 27.12 MHz,
 - Pulse frequency: 50, 100, 200, 400, 800 Hz,
 - Pulse width (duration): 20, 40, 65, 100, 200, 400 μ s,
- Foiled control panel,
- Mobile on 4 double casters including 2 with brake,
- Included:
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 45 x 52.5 x 86 cm,
- Weight: 25 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

