



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **HEINE HK 150-2 ENDO ENT Kit (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Bulb power: 15 V 150 W,
- HK 150-2 ENDO cold light projector is designed for instruments requiring fiber optic power supply. Equipped with a special halogen bulb with extended life and 150 Watt power,
- High power halogen bulb. Gives a maximum bright image.
- Aspherical lighting optics. Intense light with high color temperature,
- Infinitely adjustable light intensity without affecting the color temperature, Allows the operator to reduce the brightness of the light with the removal of red glow,
- Included:
  - Power cord,
- Dimensions: 20 x 24 x 20 cm,
- Weight: 6 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

