



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDEN-INMED AQUAI Membrane massage bed (Reconditioned)

- Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Real photos of the product,
 - Made in Poland, (2021),
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - Power: 3600 VA,
 - Aquai bed is used to perform dry hydromassage. Aquai has a wide range of applications responsible for proper physical fitness and well-being. It is perfect for both rehabilitation centers and wellness centers,
 - During the treatment, the streams of water, circulating in a closed circuit, hit the flexible membrane, so the patient can experience the healing power of hydromassage without having direct contact with the water,
 - The Aquai dry massage tub is equipped with functions for adjusting the duration of the treatment, frequency and sequence of changing massage zones. The 12 jets found in the Aquai are divided into 6 massage zones. The water temperature is also adjustable. To ensure comfort during the treatment, the water temperature is maintained throughout the treatments,
 - Basic massage programs:
 - circular massage,
 - parallel massage,
 - lymphatic massage,
 - point massage,
 - smoothing massage,
 - vibration massage,
 - pulsation massage,
 - Features:
 - Temperature selection range within: 30 °C - 40 °C,
 - 6 independent massage zones,
 - all jets can work simultaneously,
 - 6 factory massage programs,
 - possibility to set and save user programs,
 - patient comfortable piercing
 - comfortable for the patient the course of the massage with intimacy,
 - short time of preparing the patient for the procedure,
 - hygienic conditions of performing the massage,
 - does not require connection to the water and sewage system,
 - does not require drainage of water used after each treatment,
 - water heater to maintain a constant temperature,
 - touch control panel,
 - LED backlighting at the base,
 - Includes:
 - User manual in Polish (PDF),
 - Power cord,
 - Movable headrest,
 - Lever with support,
 - Hose and fittings for filling,
 - Water preservative,
 - Capacity: 220 l,
 - Dimensions:
 - Overall: 115 x 240 x 65 cm,
 - Membrane area: ~2.4 m²
 - Maximum lifting capacity: 135 KG,
 - Weight:
 - Without water filling: 190 kg,
 - With water fill: 410 kg
 - Has a current inspection and is ready for operation,
 - Issued Technical Passport (Service Report) valid for a period of 12 months,
 - Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
 - Financing options (Poland only): Installments, Leasing, Loan,
- If you have any questions, please do not hesitate to contact us!
- In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

